


# Sani-Cloth® Wipes Materials Compatibility Reference Guide



**COMPATIBLE**  
by design™




Distributed by Amtech Medical

**PDI**  
BE THE DIFFERENCE

# PDI's Compatibility Approach.

Our scientists use a variety of ASTM and ISO material testing standards to understand the effects of our disinfectants on common material surfaces.



## Wipe test:

Material surfaces were wiped 120 times with a 30-minute dry period between each wipe. Simulates time-dependent material exposure in a healthcare setting.

## Wipe test under constant strain:

Materials were wiped 120 times with a 30-minute dry period between each wipe under 1% constant strain. Simulates time-dependent material exposure under strain in the healthcare setting.



## Immersion test:

Coupons of the materials were wrapped inside the disinfectant wipe and placed inside Ziploc-type bags. Three times each week, the test samples were unwrapped and inspected for deleterious effects. Observations were noted and fresh wipes were used to wrap the test surfaces and continue the exposure. This method examines the durability of the material under extreme chemical exposure.

We have partnered with industry leaders in raw material manufacturing for medical devices to evaluate our broad family of products on thermoplastic chemical resistance.<sup>1,2,3,4,5</sup>

## PDI's Compatible by Design™ System.

### Legend



No visible surface change (damage/degradation) to the material is likely to occur when used according to the directions for use.



Some visible effects or cosmetic changes such as streaking, discoloration, or clouding, may be observed with long-term exposure.



Visible surface change (damage/degradation) may occur with long-term exposure.

# Material Compatibility Reference Guide




|                      | Surface material                      | Found in   | Super Sani-Cloth® Wipes | Sani-Cloth® Prime Wipes | Sani-Cloth® Bleach Wipes | Sani-Cloth® AF3 Wipes | Sani-Cloth® Plus Wipes | Sani-HyPerCide® Wipes |
|----------------------|---------------------------------------|--|-------------------------|-------------------------|--------------------------|-----------------------|------------------------|-----------------------|
| Polymers             | Polymethyl Methacrylate (PMMA)        | Incubators, infusion pumps, phone displays, dialyzers, X-ray protective shields  | ☺☺☺                     | ☺☺                      | ☺                        | ☺☺☺                   | ☺☺                     | ☺☺                    |
|                      | High Density Polyethylene             | Chemical containers, cable insulation, packaging, electrical device enclosures, trays, industrial plastic products   | ☺☺☺                     | ☺☺                      | ☺☺                       | ☺☺                    | ☺☺                     | ☺☺                    |
|                      | Acrylonitrile Butadiene Styrene (ABS) | Keyboards, infusion pumps, inhalers, tracheal tubes, ventilator valves, electrical device enclosures, medical masks, medical devices for blood access                                      | ☺☺☺                     | ☺☺                      | ☺☺☺                      | ☺☺                    | ☺☺                     | ☺☺                    |
|                      | Polypropylene                         | Surgical trays, device exteriors, membrane oxygenators, tubing connectors, plastic bottles   | ☺☺☺                     | ☺☺☺                     | ☺☺                       | ☺☺                    | ☺☺                     | ☺☺                    |
|                      | Polycarbonate                         | Infusion pumps, anesthesia containers, IV lock box, safety goggles, hemodialysers, blood oxygenators, blood reservoirs, surgical instruments, endoscopic appliances, lenses, IV connectors | ☺☺                      | ☺☺                      | ☺☺                       | ☺☺                    | ☺☺                     | ☺☺                    |
|                      | Marlite® <sup>6</sup>                 | Wall panels  | ☺☺☺                     | ☺☺☺                     | ☺☺☺                      | ☺☺☺                   | ☺☺☺                    | ☺☺☺                   |
|                      | Polyvinyl Chloride                    | Blood and dialysis bags, IV bags, tubing, oxygen masks, catheters, floors, electrical cords, furniture, mattress covers  | ☺☺☺                     | ☺☺☺                     | ☺☺☺                      | ☺☺☺                   | ☺☺☺                    | ☺☺☺                   |
|                      | Tritan™ <sup>7</sup> Copolyester      | Clear device components  | ☺☺☺                     | ☺☺                      | ☺☺☺                      | ☺☺☺                   | ☺☺☺                    | ☺☺☺                   |
|                      | 3Form™ Varia™ <sup>8</sup>            | Formable acrylic sheets  | ☺☺☺                     | ☺☺☺                     | ☺☺☺                      | ☺☺☺                   | ☺☺☺                    | N/A                   |
| Fabrics              | Polyurethane                          | Upholstery, hospital bedding, head supports, armrests, surgical drapes, gaskets, lights, tubing, mattress covers   | ☺☺☺                     | ☺☺☺                     | ☺                        | ☺☺☺                   | ☺☺☺                    | ☺☺☺                   |
| Metals               | Brass C260                            | Scalpels, scissors, pipe fittings  | ☺☺                      | ☺☺                      | ☺                        | ☺                     | ☺☺☺                    | ☺                     |
|                      | Aluminum 6061                         | Stethoscopes, scalpel, surgical trays, isolation carts, containers, gurneys, seating, monitors, handrails, walkers, wheel chairs, surgical lights  | ☺☺                      | ☺☺                      | ☺☺☺                      | ☺☺                    | ☺☺                     | ☺☺                    |
|                      | Galvanized Steel                      | Carts, ductwork, pipes, nails, bolts   | ☺☺☺                     | ☺☺☺                     | ☺                        | ☺☺                    | ☺☺☺                    | ☺                     |
|                      | Stainless Steel 316                   | OR tables, case carts, cabinets, carts, trolleys, wheelchairs, sinks, thermometers, hemostats, tweezers, forceps, bed frames, furniture, fixtures, counters                                | ☺☺☺                     | ☺☺☺                     | ☺                        | ☺☺                    | ☺☺                     | ☺☺                    |
|                      | Aluminum Silicate                     | Fixtures, insulators, Instrument trays, walkers, crutches, folding stretchers  | ☺☺☺                     | ☺☺☺                     | ☺☺                       | ☺☺☺                   | ☺☺☺                    | ☺☺☺                   |
| Hard Porous Surfaces | Ceramic                               | Pressure sensors, hand tools, valves, tiles  | ☺☺☺                     | ☺☺☺                     | ☺☺☺                      | ☺☺☺                   | ☺☺☺                    | ☺☺☺                   |
|                      | Porcelain                             | Bathrooms  | ☺☺                      | ☺☺                      | ☺☺☺                      | ☺☺☺                   | ☺☺                     | ☺☺☺                   |
|                      | Marble                                | Decorative countertops   | ☺☺☺                     | ☺☺☺                     | ☺☺☺                      | ☺☺☺                   | ☺☺                     | ☺☺☺                   |

## The challenge: equipment damage.

Incompatible cleaning and disinfection products in healthcare creates dangerous safety issues for patients, costs healthcare facilities millions of dollars annually, and impedes effective cleaning and disinfection.


Biomedical engineers estimate a 20-27% equipment life span reduction due to incompatible disinfectants.<sup>9</sup>



## Types of surface change.

Exposure to disinfectants can lead to degradation of surfaces and, for plastics, environmental stress cracking (ESC).

### Discoloration



Material's surface changes color due to environmental conditions. Often described as fading, yellowing, staining, darkening, bleaching, spotting, etc.

### Residue



Material remaining on the surface after the solvent or additives have evaporated. Residue is unsightly but often removable.

### Cracking




Material failure, fracture, crack, and crazing described as localized deformation occurring in a material.

### Corrosion



Gradual deterioration of materials caused by a chemical reaction on the surface. Can lead to pitting, rusting, deformation, etc.

### Delamination



Mode of failure where a material fractures into layers along a plane parallel to a surface.

## Guidance from the CDC, FDA and Joint Commission.

“Disinfect noncritical surfaces and medical devices with an EPA-registered hospital disinfectant using the label’s safety precautions and use directions.”<sup>10</sup>

- CDC from “Guideline for Disinfection in Healthcare facilities”

“FDA recommends that you validate your disinfection processes and instructions. FDA also recommends that you follow the recommendations in device-specific FDA guidance documents or any relevant FDA-recognized standards.”<sup>11</sup>

- FDA from “Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling”



The goal of our Compatible by Design Program is to help your facility meet the requirements of your accreditation organization.





# Common surface materials found in a patient room.



## PDI's Compatibility Search Tool.

WE ARE PDI

UNITED STATES ABOUT PDI CAREERS

HOME VIEW PRODUCTS WHO WE SERVE PDI EDUCATION CUSTOMER CARE CONTACT US

### Healthcare Professionals.

Do you know your medical devices' disinfection guidelines? We can help identify which products are approved and compatible!

[CLICK HERE TO RECEIVE MORE INFORMATION](#)

#1 BRAND OF DISINFECTANT WIPES IN HEALTHCARE

Check your Equipment Compatibility.  
Select a PDI Product.

Manufacturer

Equipment Type

PDI Product



Scan here to learn more about PDI's Compatibility Search Tool or go to <https://pdihc.com/compatibility-hcp>

# PDI's compatible disinfectants.

## Sani-Cloth® Prime Germicidal Disposable Wipe



Powered by a next generation formulation effective against HAI causing pathogens in just 1 minute.



## Sani-HyPerCide® Germicidal Disposable Wipe



A powerful, ready to use hydrogen peroxide formula designed to address HAI causing microorganisms, including *Clostridioides difficile*.



## Super Sani-Cloth® Germicidal Disposable Wipe



Ideal for daily use in fast-paced environments that require a short contact time and broad coverage of microorganisms.



## Sani-Cloth® AF3 Germicidal Disposable Wipe



Alcohol free. Fragrance free. Worry free. Ideal for use around patients and staff, especially those with respiratory sensitivities.



## Sani-Cloth® Plus Germicidal Disposable Cloth



A low-level disinfecting wipe for general surface disinfection.



## Sani-Cloth® Bleach Germicidal Disposable Wipe



Ideal for disinfecting areas contaminated with *Clostridioides difficile* spores and Norovirus.



# Educational and point-of-care accessories.

Developed to assist you in facility protocol compliance.



PDI Sani-Bracket® for wall/pole mounting of a canister



Sani-Canister Caddy® for placement on flat surfaces (bilingual)



Sani-TAG® Equipment ID System for identification of equipment to ensure the proper wipe is used on a particular piece of equipment



Educational Wall Charts and signage (bilingual)



# Order information.

REORDER NO. WIPE SIZE CASE PACK CASE WGT CASE CUBE PALLET TI/HI



## Sani-Cloth® Prime Germicidal Disposable Wipe

Powered by a next generation formulation, help protect your patients, staff, and facility from HAIs in just 1 minute.

|                           |        |                |         |           |          |      |
|---------------------------|--------|----------------|---------|-----------|----------|------|
| Large Canister            | P25372 | 6" x 6.75"     | 12/160s | 30.47 lbs | 1.430 ft | 10/3 |
| X-Large Canister          | P24284 | 7.5" x 15"     | 6/70s   | 18.44 lbs | 0.936 ft | 10/4 |
| Large Individual Packet   | H06182 | 5" x 8"        | 10/50s  | 9.99 lbs  | 0.636 ft | 17/4 |
| X-Large Individual Packet | U13195 | 11.5" x 11.75" | 3/50s   | 8.94 lbs  | 0.500 ft | 24/3 |



## Sani-HyperCide® Germicidal Disposable Wipe

A ready to use hydrogen peroxide formula designed to address HAI causing microorganisms, including *Clostridioides difficile* without compromising compatibility.

|                  |        |            |         |           |          |      |
|------------------|--------|------------|---------|-----------|----------|------|
| Large Canister   | P27372 | 6" x 6.75" | 12/160s | 26.49 lbs | 1.430 ft | 10/3 |
| X-Large Canister | P26584 | 7.5" x 15" | 6/65s   | 15.04 lbs | 0.936 ft | 10/4 |



## Super Sani-Cloth® Germicidal Disposable Wipe

2 minute contact time allows for a quick room turnover. Ideal for daily use in fast-paced environments that require short contact times and broad coverage of microorganisms.

|                           |        |                |         |           |          |      |
|---------------------------|--------|----------------|---------|-----------|----------|------|
| Large Canister            | Q55172 | 6" x 6.75"     | 12/160s | 26.04 lbs | 1.430 ft | 10/3 |
| X-Large Canister          | P86984 | 7.5" x 15"     | 6/75s   | 15.56 lbs | 0.936 ft | 10/4 |
| Large Individual Packet   | H04082 | 5" x 8"        | 10/50s  | 8.03 lbs  | 0.539 ft | 20/4 |
| X-Large Individual Packet | U87295 | 11.5" x 11.75" | 3/50s   | 7.19 lbs  | 0.500 ft | 24/3 |
| Softpack                  | A22480 | 8.2" x 9.8"    | 9/80s   | 16.58 lbs | 0.948 ft | 10/4 |



## Sani-Cloth® AF3 Germicidal Disposable Wipe

Helping to protect your staff, patients, and residents has never been easier! Alcohol and fragrance-free formula is ideal for use around those with respiratory sensitivities.

|                           |        |                |         |           |          |      |
|---------------------------|--------|----------------|---------|-----------|----------|------|
| Large Canister            | P13872 | 6" x 6.75"     | 12/160s | 26.04 lbs | 1.430 ft | 10/3 |
| X-Large Canister          | P72584 | 7.5" x 15"     | 6/75s   | 16.71 lbs | 0.936 ft | 10/4 |
| Large Individual Packet   | H59200 | 5" x 8"        | 10/50s  | 8.14 lbs  | 0.539 ft | 20/4 |
| X-Large Individual Packet | U27500 | 11.5" x 11.75" | 3/50s   | 7.67 lbs  | 0.500 ft | 24/3 |
| Pail                      | P1450P | 7.5" x 15"     | 2/160s  | 12.40 lbs | 0.953 ft | 10/5 |
| Refill (for Pail)         | P2450P | 7.5" x 15"     | 2/160s  | 11.10 lbs | 0.529 ft | 18/4 |



## Sani-Cloth® Plus Germicidal Disposable Wipe

Ideal for use in alternate care settings, including physician and dental offices.

|                  |        |            |         |           |          |      |
|------------------|--------|------------|---------|-----------|----------|------|
| Large Canister   | Q89072 | 6" x 6.75" | 12/160s | 26.14 lbs | 1.430 ft | 10/3 |
| X-Large Canister | Q85084 | 7.5" x 15" | 6/65s   | 14.85 lbs | 0.936 ft | 10/4 |



## Sani-Cloth® Bleach Germicidal Disposable Wipe

Ideal for disinfecting areas contaminated with *Clostridioides difficile* spores and Norovirus.

|                              |        |                |         |           |          |      |
|------------------------------|--------|----------------|---------|-----------|----------|------|
| Clinical Size Wipes Canister | P84172 | 6" x 5"        | 12/160s | 26.86 lbs | 1.430 ft | 10/3 |
| Large Canister               | P54072 | 6" x 10.5"     | 12/75s  | 26.50 lbs | 1.430 ft | 10/3 |
| X-Large Canister             | P25784 | 7.5" x 15"     | 6/65s   | 19.79 lbs | 0.936 ft | 10/4 |
| Large Individual Packet      | H58195 | 5" x 7"        | 10/40s  | 8.25 lbs  | 0.505 ft | 19/5 |
| X-Large Individual Packet    | U26595 | 11.5" x 11.75" | 3/40s   | 8.20 lbs  | 0.500 ft | 24/3 |
| Pail                         | P7007P | 7.5" x 15"     | 2/160s  | 16.50 lbs | 0.953 ft | 10/5 |
| Refill (for Pail)            | P700RF | 7.5" x 15"     | 2/160s  | 15.60 lbs | 0.529 ft | 18/3 |



Scan here to learn more about PDI's compatibility partnerships or go to <https://pdihc.com/equipment-compatibility>

<sup>1</sup>Data on file. Compatibility study 8-23.2017.

<sup>2</sup>Zettel, M., Liu, Y., Lamont, M., Bernhard, K. and Hicks J., October 2022. "Standard test method for evaluation of chemical resistance of thermoplastics for non-disposable medical devices- a voice from industry collaboration through the COVID-19 pandemic". Medical Design Briefs. <https://www.medicaldesignbriefs.com/component/content/article/mbd/pub/features/articles/46685>.

<sup>3</sup>Esposito, A. October 2016. "PDI and SABIC join forces to test medical devices to improve patient protection from infection". Medical Design & Outsourcing. <https://www.medicaldesignandoutsourcing.com/pdi-sabic-join-forces-test-medical-device-plastics-help-improve-patient-protection-infection/>

<sup>4</sup>Ratnagiri, R., and Martin, R. (2022). Resistance of DuPont polymers to disinfecting chemicals, including those effective against the SARS-CoV-2 virus. [White paper], Dupont. <https://www.dupont.com/content/dam/dupont/amer/us/en/mobility/public/documents/en/Home-Appliance-Whitepaper-FNL.pdf>

<sup>5</sup>(2020). Covestro Polycarbonates: Compatibility with Disinfectants used against SARS-CoV-2. [White paper], Covestro. <https://solutions.covestro.com/-/media/covestro/solution-center/whitepapers/cov-chemical-compatibility-to-disinfectants-used-against-sars-cov-2-2020-06.pdf>

<sup>6</sup>Marlite® is a registered trademark of Marlite, Inc.

<sup>7</sup>Tritan™ is a trademark of Eastman Chemical Company.

<sup>8</sup>3Form™ and Varia™ are trademarks of 3Form, LLC

<sup>9</sup>Data on file. Compatibility Study 8-28-17.

<sup>10</sup><https://www.cdc.gov/infectioncontrol/guidelines/disinfection/index.html#rec5g>

<sup>11</sup><https://www.fda.gov/regulatory-information/search-fda-guidance-documents/reprocessing-medical-devices-health-care-settings-validation-methods-and-labeling>

<sup>12</sup>Sani-Cloth® Plus Germicidal Disposable Cloth compares to competitor quat/low-alcohol disinfectant products.

©2023 PDI 400 Chestnut Ridge Road Woodcliff Lake, NJ 07677 T: 800-999-6423 W: [pdihc.com](https://pdihc.com)  
PDI09238824



Distributed by  **amtech**

Care without compromise.

0800 268 324 | [amtech.co.nz](https://amtech.co.nz)

**Auckland Office** 525 Great South Road, Penrose, Auckland  
**Head Office** 281 Heads Road, PO Box 2059, Wanganui